

---

# A Comparison of Methods for the Detection of Smooth Surface Caries

B. Wong<sup>1</sup>, K. Sivagurunathan<sup>1</sup>, J. D. Silvertown<sup>1</sup>, W. M. P. Hellen<sup>2</sup>, G. I. Elman<sup>2</sup>,  
**S. H. Abrams<sup>1</sup>**, L. O. Okoye<sup>3</sup>, B. T. Amaechi<sup>4</sup>

<sup>1</sup>Quantum Dental Technologies, Toronto, ON, Canada

<sup>2</sup>Cliffcrest Dental Office, Scarborough, ON, Canada

<sup>3</sup>University of Nigeria, Enugu, Nigeria

<sup>4</sup>Comprehensive Dentistry, University of Texas Health Science Center at San Antonio, San Antonio, TX


# Introduction

---

- Detection of non-cavitated caries is important because lesion progression may be halted at this stage, remineralized or minimally restored, thereby preserving natural tooth structure.
- Visual and tactile methods of caries detection only examine the tooth surface and not the lesion developing beneath it.
- Smooth surface changes may be detected visually but are there other methods to detect & monitor lesion changes over time?

# Objective

This *in vitro* study evaluated the ability of The Canary System, DIAGNOdent, Spectra Caries Detection Aid, ICDAS II, and Radiographic Examination to detect smooth surface caries.



**Spectra Caries  
Detection Aid**



**The Canary System™**



**DIAGNOdent**



**Radiographic Examination**

#### Caries Codes

- 0 = Sound tooth surface
- 1 = First visual change in enamel
- 2 = Distinct visual change in enamel
- 3 = Enamel breakdown, no dentine visible
- 4 = Underlying dentinal shadow (not cavitated into dentine)
- 5 = Distinct cavity with visible dentine
- 6 = Extensive distinct cavity with visible dentine

**Visual Examination (ICDAS II)**

# Materials and Methods

---


- 92 healthy & carious sites on smooth surfaces of extracted human teeth were used.
- A blinded, experienced operator scanned the teeth using The Canary System, DIAGNOdent and Spectra with three repeat measurements per site.
- Two blinded clinicians independently scored the teeth using ICDAS II.
- The same two blinded clinicians independently ranked radiographs of the teeth as '1' for presence of caries and '2' for absence of caries. **Note radiographs were taken with smooth surfaces mounted as interproximal lesions**
- Where there was disagreement between the clinicians' scores, the site were re-examined by both clinicians together and a consensus score reached.

# Canary Scale

---




Source: The Canary System User Manual




# Examining Lesions with Canary

---



- Angulation of Canary Tip will provide a range of Canary Numbers depending upon what is beneath the beam
- As one scans along the occlusal surface one can detect & image the lesion.
- The Canary acts like a “punch biopsy” for examining the tooth surface



CANARY  
NUMBER


16

**Scanning the Occlusal Surface to Map the Lesion**






CANARY NUMBER
16
19




CANARY NUMBER
16
19
35

**Scanning the Occlusal Surface to Map the Lesion**



CANARY NUMBER
16
19
35
75

**Scanning the Occlusal Surface to Map the Lesion**



CANARY NUMBER
16
19
35
75
80

**Scanning the Occlusal Surface to Map the Lesion**

# DIAGNOdent Scale

---

0-10	Healthy Tooth Structure
11-20	Outer Half Enamel Caries
21-30	Inner Half Enamel Caries
30+	Dentin Caries

Source: DIAGNOdent Operating Guide



# Materials and Methods - Validation

---

- Polarized Light Microscopy (PLM) was performed blinded at the Department of Comprehensive Dentistry, University of Texas Health Science Center at San Antonio as validation.

# Statistical Analysis - Correlation

---

- Correlation between ICDAS II scores and the numerical readings from The Canary System, DIAGNOdent and Spectra and the scores from Radiographic Examination were determined by Pearson's coefficient of correlation ( $R^2$ ,  $p < 0.01$ ).
- Correlation between lesion depth and the numerical readings from The Canary System, DIAGNOdent, Spectra and ICDAS II scores were determined by Pearson's coefficient of correlation ( $R^2$ ,  $p < 0.01$ ).



# Statistical Analysis – Sensitivity & Specificity

---

Sensitivity and specificity were determined using PLM results and the following criteria:

Device	Sound	Caries
Canary Number	$\leq 20$	$> 20$
DIAGNODent	$\leq 10$	$> 10$
SPECTRA	$\leq 1$	$> 1$
ICDAS II	$= 0$	$\geq 1$

---

# RESULTS

# Correlation with ICDAS II

---

Caries Detection Method	Correlation with ICDAS II Scores (R <sup>2</sup> )
The Canary System	0.798
DIAGNOdent	0.244
Spectra	0.592
Radiographic Examination	0.091

# Correlation with Lesion Depth (PLM)


---

Caries Detection Method	Correlation with Lesion Depths ( $R^2$ )
The Canary System	0.583
DIAGNOdent	0.550
Spectra	0.423
ICDAS II	0.470

# Representative Sample (#5A)


Representative sample with visually-carious examination site A and healthy examination site B.

Site A  
ICDAS  $\geq 1$




True Positive

Site A  
Canary Number = 43




True Positive

Site A  
DIAGNOdent = 57




True Positive

Site A  
Spectra Value = 0




False Negative

Site A  
PLM = 874  $\mu\text{m}$   
Gold Standard






# Radiographic Exam Sample 5A

---




# Sample 10B



Device	Reading	
Canary	66	True Positive
DIAGNODent	1	False Negative
ICDAS II	2	True Positive
SPECTRA	0.9	False Negative
PLM Depth	550.91 microns	

# Sample 10 Radiograph

---





# Sensitivity & Specificity

---

Caries Detection Method	Sensitivity	Specificity
The Canary System	0.84	0.91
DIAGNOdent	0.49	1.00
Spectra	0.51	1.00
ICDAS II	0.83	1.00

# Conclusions: Detecting Smooth Surface Caries

---

- Strong positive correlation between Canary Numbers & ICDAS II scores for detecting smooth surface caries.
- Spectra, DIAGNOdent & Radiographic Examination demonstrated poorer correlation with ICDAS II.
- ICDAS II may not be as sensitive to changes in lesion size within each classification.
- The strong correlation between The Canary System and ICDAS II implies that these two methods may be combined to increase their effectiveness for detection of caries on **smooth surfaces**.
- ICDAS II & The Canary System showed superior sensitivity compared to DIAGNOdent and Spectra.



# Thank You

For more information contact:

**Stephen Abrams, DDS**  
President, Quantum Dental Technologies

t. (416) 265-1400 (dental practice)  
m. (416) 523-8453  
e. [stephen@thecanarysystem.com](mailto:stephen@thecanarysystem.com)